Japet^W.

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Table of contents

01.	Symbols	29
	Properties and copyright	29
02.	Presentation and safety	30
	Device overview	30
	Clinical performances	30
	Indications	31
	Contra-indications	31
	Warnings	31
	Possible adverse effects of the device	32
	Vigilance	32
03.	Precautions for use	32
	Body measurements	32
	Buttons and LED description	34
	Device's operation	34
	a. Setting up the device	34
	b. Switching On	37
	c. Decompression level adjustment	37
	d. Switching Off	38
	Battery	38
04.	External presentation	39
	Equipment included	39
	Content	39
	General front view	39
	Right side view	39
	Battery access	39
05.	Maintenance and cleaning	40
	Battery charge	40
	Battery disposal	40
	Device cleaning	41

06.	Technical characteristics	41
	Device features	41
	Battery electrical features	41
	Charger electrical features	42
	Mechanical properties	42
	Claimed performance	42
	Connectivity features	42
07.	Conditions of use, storage and transport	42
	Operating conditions	42
	Storage and transport conditions	42
	Electromagnetic compatibility	43
08.	Troubleshooting	44
09.	Regulatory information	44
	Responsibility	44
	Guarantees	44
	Registered trademarks	45
	Patents	45
10.	Other information	46
	Symbols on the label and packaging	46
	Version history of the user guide	49
	Contact information	

About this guide

This user guide has no contractual value and JAPET MEDICAL DEVICES® cannot be held liable under any circumstances on the basis of information contained in this guide.

This user guide details all the necessary knowledge for the implementation, use and maintenance of the JAPET.W device and the collection of information displayed.

After a careful and complete reading before using the device, the operator will be able to:

- Connect the peripheral elements (power cord) and turn on the JAPET.W device.
- Configure the device, perform basic maintenance.

JAPET MEDICAL DEVICES® studies this guide "as is", without warranty of any kind, either express or implied, including but not limited to, the implied warranties of sales terms and/or suitability for a particular purpose in order to provide simple and accurate information. The company JAPET MEDICAL DEVICES® can therefore assume no responsibility for any or misinterpretation. Although every effort has been made to provide a manual as accurate as possible, it may contain technical inaccuracies and/or typographical errors. In no event should JAPET MEDICAL DEVICES® be liable for any loss of profit, loss of business, loss of data, business interruption, or for indirect, special, incidental, or consequential damages of any kind.

In the event of damage arising from a defect (imperfection) or error contained in this user guide, JAPET MEDICAL DEVICES® undertakes to send as soon as possible to the client a paper or electronic document containing the corrigenda mad to this

guide. This guide is regularly updated. The most recent version of this guide is available upon request from JAPET MEDICAL DEVICES®. However, in the event of major modifications being made to the guide, JAPET MEDICAL DEVICES® undertakes to send the new paper or electronic version of the guide to the client as soon as possible. It is specified that this does not imply updating the hardware and/or the software in your possession.

The owner of the product has to keep this manual for the duration of the use of the product. This guide includes a chapter on quick resolution of the most common problems.

Any request for information or modification relating to this guide should be addressed to :

JAPET MEDICAL DEVICES® ,
F3A
12, avenue Pierre Mauroy
59120 LOOS
FRANCE

01. Symbols

Warning

Attention at this level indicates potential property damage that may produce inaccurate data or prevent a device function, even if personal injury is unlikely.

FR

EN

DE

Properties and copyright

All manuals and documentation of any kind are the property of JAPET MEDICAL DEVICES® and are protected by copyright, all right reserved. Your right to copy this documentation is limited to legal copying rights. These manuals may not be distributed, translated or reproduced, in whole or in part, in any manner or form, without the prior written permission of JAPET

MEDICAL DEVICES®. Reproduction, adaptation or translation of this guide without prior written permission is prohibited, to the extent permitted by copyright laws. Copyright® – 06/2018 – JAPET MEDICAL DEVICES® – All rights reserved.


02. Presentation and safety

Device overview

Low back pain is a major and growing public health problem in Western countries where it affects or will affect nearly 4 out of 5 people. Among its main causes, the squashing of intervertebral discs can cause from a simple persistent pain to paralysis. This pain, when it persists and intensifies, creates a disabling pain that restricts the physical activities and progressively isolates from the social and professional environment. The complexity of the origin of the pathology and of the multidisciplinary management amplifies the risks of transition to chronicity.

JAPET.W, the active exoskeleton to relieve pain related to low back pain, while allowing a complete mobility.

It is inspired by devices commonly used in rheumatology such as elongation tables to relieve the patient by decompressing the spine. The traction is produced by the thrust of micromotors in between on two belts, one on the pelvis, the other on the upper body.

The device is also inspired by robotic exoskeletons in order to adjust the position of the micromotors according to the user movements to release trunk mobility. He can thus resume a recommended physical activity for this type of pathology.

Clinical performance of the device.

The Japet.W provide a decompressions in the intervertebral disc.

Indications

The JAPET.W device is indicated for people suffering from chronic and acute low back pain. The device is intended for ambulatory distension of the lumbar spine.

It is recommended to use the device above a T-shirt.

The use of the device is intended for medical professionals and/or referent persons who are trained by JAPET MEDICAL DEVICES® or by an authorised third to realise this training.

It is recommended to use the device on sessions of 30 to 60 minutes (2 to 3 times a day).

For a use in the work environment, it is recommended to wear the JAPET.W device before a painful event.

Warning

Do not use batteries other than those distributed by JAPET MEDICAL DEVICES®

Do not use any battery charger other than the one distributed by JAPET MEDICAL DEVICES®

Do not use JAPET.W device if it is damaged in any way.

Take common precautionary measures to avoid any contact with blood or other biohazardous materials. Contaminated material must be treated in accordance with the hygiene requirements in force in the establishment.

Do not expose this device to rain, dripping or splashing. Do not immerse this appliance.

Contra-Indications

This device is contraindicated in the following cases:

- Medullary and radicular neurological deficits
- Neurostimulation implant therapy
- Cardiac or circulatory disease and severe respiratory problems
- Tumoral or infectious process of the vertebrae
- Osteoporosis (< -2.5 ds)
- History of arthrodesis surgery
- Herniated disc surgery (< 3 months)
- Dorso-lumbar spine fracture (< 3 months)
- Recent rib fracture (< 3 months)
- Skin lesion, contusions and stretching injuries of the trunk
- Pregnancy (> 4 months)
- History of lumbar prosthesis surgery

If any of the above conditions apply to you, consult your doctor before starting to use the appliance.

The device is internally powered by a battery that is charged outside the device using the supplied charger.

It is strictly forbidden for a third party other than a technician authorized by JAPET MEDICAL DEVICES® to open or modify the device.

It is forbidden to disassemble or modify the device (this act interrupts all guarantees).

For all cleaning operations, switch off the units and remove the battery to avoid short circuits.

The 4 motors may become hot during normal use of the device, the maximum temperature that can be in contact with the skin is 105.8 °F (41°C).

Possible adverse effects of the device

Wearing the device may cause :

- Haematomas
- Irritations

In case of an adverse effect not listed and related to the device, please contact JAPET MEDICAL DEVICES ®. Contact informations are available at the end of the user guide.

Vigilance

In the event of a serious incident, i.e., any incident that directly or indirectly resulted in, or was likely to result in ;


- The death of a patient, user or any other person;
 - A serious temporary or permanent deterioration in the state of health of a patient, user or any other person; s
 - A serious threat to public health;
- Please contact JAPET MEDICAL DEVICES ®, contact information is available at the end of the user guide


03. Precautions for use

Body measurements

There are two different device sizes, rated T1 for the smallest and T2 for the largest. In order to select the correct device size, measure the under breast circumference as well as the hip circumference as shown below:

Now enter the values in the following graphic to find the size of the device adapted to the user's morphology:

Buttons and LED description

	ON / OFF	Turning on the device / Turning off the device (an audio signal confirms the pressing)
	INCREASE	Increase decompression force (an audible signal confirms the pressure)
	DECREASE	Decrease decompression force (an audible signal confirms the pressure)
	STRENGTH INDICATOR	Indicates the current level of decompression (4 levels) Flashing: indicates an error or an empty battery.

Device's operation

Follow these steps consecutively to use the JAPET.W device.

a. Setting up the device

To position the device correctly on the user, follow the steps.

Positioning of the belt (fig.a)

1. Hold the device on the blue cases then pass it behind your back, get help if necessary.
2. The arrows indicate the lower part to position on the hips.
3. Close the lower velcro® so that the device does not fall out during adjustment, then close the upper velcro.

Tighten the belt (fig.b)

1. Install the device symmetrically on the hips. Check the alignment of the lines on the high upper and lower belts.

By pressing the arrows on the outside should feel symmetrical pressure on the hips.

2. Check that you can pass a finger between the top and blue boxes.

3. To tighten the device, press the reel (until you hear a click). Once pushed in, turn it in a clockwise direction. Perform this manipulation for the upper and lower reels.

Belt adjustment (fig.c)

1. Adjust the tightness with the reels: push the reel in until you hear a "click" and then turn it clockwise until it is difficult for you to put a finger between your body and the device, see the picture.


In order to have a real funnel effect and a good support, it is necessary to tighten the device well.

2. To loosen the device, pull the reel until you hear a "click".


The white lines on the side of the device must be aligned.

3. Do not hesitate to repeat the adjustment several times to ensure that the device is symmetrical and correctly tightened.

1


**Positioning of the belt
(fig.a)**

2

3

2

3


**Tighten of the belt
(fig.b)**

1

**Belt adjustment
(fig.c)**

2

b. Switching on

Press the “  ” for about one second

Following this, two witnesses are activated:

- An indicator light indicates the battery level temporarily. The gauge lights up gradually on all four LEDs, see table on page 19.
- An increasingly high-pitched sound signal.

c. Decompression level adjustment


Four decompression levels are available to adapt to the different needs of the users and healthcare professionals, depending on the use of the device and the exercises performed by the user.

Press the “+” or “-” button to adjust the decompression level. An audio signal confirms the procedure.

Refer to the table below:

 If the device is not used over a long period of time, the battery must be taken out of the device.

 If the decompression indicator flashes, it is an operating error.

	Level 4	Maximum decompression
	Level 3	Upper intermediate decompression
	Level 2	Lower intermediate Decompression
	Level 1	Normal decompression

d. Switching off

1. First press the “  ” button up to threshold 0, the actuators (blue cases) retract halfway, or completely if the actuators are slightly extended.

2. The device can be removed.

3. Press the “  ” button for one second to deactivate the device.

4. The actuators will return to the initial position and all indicators will turn off.

When the battery is inserted, the actuators retract automatically if they are not in their initial position.

Recharge the battery by following the instructions in the maintenance and cleaning>recharging section.

In case of persisting problems, contact JAPET MEDICAL DEVICES®.

Battery

An indicator light indicates the battery level.

The gauge lights up gradually on the four white LEDs for 1 second.

Refer to the table below:

Contact information is available at the end of the user guide.

	Full battery	Four LEDs light up gradually, battery level is 100%.
	Batterie almost full	Three LEDs light up gradually, battery level is 75%
	Battery half depleted	Two LEDs light up gradually, battery level is 50%.
	Battery almost depleted	Only one LED lights up, the battery level is 25%.
	Battery depleted	Only one LED flashes, battery level is below 5%. Recharge the battery according to the instructions.

04. External presentation

FR

EN

DE

Equipment included


If the protective case is damaged, please contact JAPET MEDICAL DEVICES®. When unpacking the package, check the conformity of the contents with the following list:

- 1 JAPET.W device
- 1 User guide
- 2 Batteries - RRC 2040
- 1 Charger and its power cord


Warning : To avoid any electrical short-cut risk, the charger must only be connected with the supplied base socket with an earth connection for protection.

Content

General front view

Right side view

Battery acces


05. Maintenance and cleaning

Battery charge

Battery charge When the battery level is insufficient (about 5% remaining), two indicators are activated:


- Persistent high-pitched beep
- One indicator light (flashing) at the first white LED. (fig.a)

After this, remove the battery from the device and place it in the charger as shown below. (fig.b)

Please ensure that the charger is connected to the mains supply using the supplied cable.

- A red light indicates that the battery is charging.
- A green light indicates that the charge is complete.

Do not use batteries other than those distributed JAPET MEDICAL DEVICES®

Contact information is available at the end of the user guide.

Battery disposal

Directive 2006/66/EC of the European Parliament and of the Council of 6 September 2006 on batteries and accumulators and waste batteries and accumulators regulates the recycling of batteries, electrical wires and electronic components.

The JAPET.W battery contains hazardous products, such as lithium, which cannot be disposed of with household waste. Electrical wires and other electronic components may not be disposed of with household waste.

Do not dispose of the batteries supplied by JAPET MEDICAL DEVICES® with normal household waste at the end of the product's life, but dispose of them in an appropriate recycling centre, therefore helping to protect the environment.

Device cleaning

After each use, JAPET MEDICAL DEVICES[®] recommends the use of a fabric complaint bactericidal and yeast spray, for cleaning the

device. For better results, gently wipe with a cloth.

Do not clean the device when the battery is connected.

Do not clean the device with a washing machine or by hand.

06. Technical characteristics

Technical specification

Manufacturer:	JAPET MEDICAL DEVICES [®] FA3 – 12, avenue Pierre Mauroy 59120 LOOS FRANCE
Model:	JAPET. ^W
Medical device classification:	Class IIa according to Rule 9, Annex IX of Directive 93/42/EEC
Type of applied part	BF (whole JAPET.W device)
IP Code	IP21
Embedded software class	Classe A
Service life	5 years

Battery electrical features

Model	RRC2040
Rated voltage	10,8V
Capacity	3350mAh
Dimensions	84,9mm x 58,8mm x 21,9mm
Weight	170 g
Lifetime at 77°F	> 300 cycles with a minimum of 75% of its initial capacity

Charger electrical features

Model	RRC-SMB-MBC
Input voltage range	100 - 240VAC
Input current	2,80A max
Stand by power	No load <0.5 W @ 230 VAC

Mechanical properties

Dimensions	100cm x 30cm x 10cm
Weight of the device	< 2 Kg (without battery)

Claimed performance

Force per actuator	4Kg
Global force	16Kg

Connectivity features

Connectivity	Bluetooth BLE
Frequency bandwidth (receiving)	2400 MHz à 2483,5 MHz
Frequency bandwidth (emitting)	2400 MHz à 2483,5 MHz
Type	Modulation BLE
Radiated power	< 0 dBm

07. Conditions of use & storage and transport

Operating conditions

Temperature:	+32°F à +104 °F (+0°C à +40°C)
Humidity in the air:	5% to 95%

Storage and transport conditions

Temperature:	-4°F to +104°F (-20°C to +40°C)
Humidity in the air:	5% to 95%
Atmospheric pressure:	57 kPa to 106 kPa

Electromagnetic compatibility

Electromedical devices are subject to special electro-magnetic compatibility (EMC) precautions and should only be installed and used in accordance with the information specified in this document. Wireless communication devices such as wireless home network routers, mobile phones, cordless phones and their bases and walkie-talkies may interfere with JAPET.W. Other cables and accessories not supplied by the manufacturer may affect EMC performance and cause increased emissions or decreased electromagnetic immunity of JAPET.W.

The JAPET.W device must be set up and put into use in hospitals away from high frequency MRI imaging source.

The performances of the JAPET.W device amount an overall traction force of up to 16 kg. In case of electromagnetic perturbations, these performances might not be reached and the therapeutic effectiveness of the JAPET.W limited.

The JAPET.W device is compliant with the IEC 60601-1-2: 2014 (Ed. 4).

Warning: it is recommended to avoid using this device right next to other devices or stacked on top of another since that could lead to potential malfunction. In case that would be necessary, it is recommended to monitor this device and the others to assess their proper functioning.

Warning: The use of accessories, transducers or cables other than those specified or provided with this device can lead to an increase of the electromagnetic emissions or a decrease of the immunity of this device and result in an improper functioning.

Warning: It is recommended not to use any RF portable communication device (including devices such as antenna cables or external antennas) less than 12 inches from any part of the JAPET.W device. Otherwise, the performances of this device could be altered.

08. Troubleshooting

If an error occurs, an audible signal is emitted continuously and the decompression indicator flashes.

Press and hold the “” button on the device to turn it off.

Press the “” button again to restart the device.

If the error persists, remove the battery, make sure it is charged, then return it to the device.

Then press the “” button to restart the device.

In case of persistent problems, contact JAPET MEDICAL DEVICES®.

Do not disassemble the device.

In case of requesting a new battery, contact JAPET MEDICAL DEVICES®.

09. Regulatory informations

Responsability

JAPET MEDICAL DEVICES® cannot be held responsible for any damage that may result from the use of this product under conditions that do not comply with the instructions given in this guide.

Garantees

JAPET MEDICAL DEVICES® warrants to the purchaser that the device («the warranted product») is free from defects in workmanship or materials under normal, appropriate and intended use.

JAPET MEDICAL DEVICES®'s obligations under this warranty are to repair or replace all or part of the warranted product that JAPET MEDICAL DEVICES® deems reasonable to cover because defective in parts and labour, provided that the purchaser asserts its warranty during the warranty period

and the product is returned to JAPET MEDICAL DEVICES® in accordance with the provisions of the general conditions of sale or distribution contracts. Repair or replacement of products under this warranty does not extend the warranty period.

To request repair or replacement under this warranty, the purchaser must contact Japet Medical Devices® directly. Japet Medical Devices® will authorise the purchaser to return all or part of the product guaranteed.

JAPET MEDICAL DEVICES® will decide whether to repair or replace the product and parts covered by this warranty and all replaced products or parts will become the property of Japet Medical Devices®.

During the warranty, Japet Medical Devices® may, but is not obligated to make engineering improvements to all or part of the warranted product. If Japet Medical Devices® deems it reasonable that a repair or replacement is covered by the warranty, Japet Medical Devices® will assume the cost of shipping the repaired or replacement product to the purchaser. Under this warranty, the risk of loss or damage during shipment shall be borne by the party shipping the product. Products returned by the purchaser under this warranty will be packaged in the original packaging or equivalent packaging to protect the product. If the purchaser returns a product to Japet Medical Devices® in an inadequate packaging, any damage to the material noted upon receipt by Japet Medical Devices® (and not already reported) will be presumed to have occurred during transit and will be the responsibility of the purchaser.

This warranty does not extend to any or all of the products covered by this warranty that have been misused, neglected, or damaged by causes external to the products, used against Japet Medical Devices® instructions, modified, disassembled, repaired, or reassembled by a technician who is not a Japet Medical Devices® employee. Japet Medical Devices® shall not be liable for any repairs, replacements or modifications resulting in whole or in part from normal wear and tear of the equipment. Japet Medical Devices® does not grant its warranty:

- (a) any product that is not warranted,
- (b) any product purchased from a supplier other than Japet Medical Devices® or its authorized distributor; or
- (c) any product sold under a brand name other than Japet Medical Devices®.

THIS WARRANTY IS THE ONLY AND EXCLUSIVE WARRANTY OF JAPET MEDICAL DEVICES® PRODUCTS GRANTED SOLELY TO THE PURCHASER AND EXPRESSLY REPLACES ANY OTHER EXPRESS OR IMPLIED WARRANTY INCLUDING, BUT NOT LIMITED TO, ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. JAPET MEDICAL DEVICES®'S LIABILITY ARISING FROM THE SALE OR USE OF THE PRODUCTS UNDER WARRANTY, CONTRACT, TORT OR OTHERWISE SHALL NOT EXCEED THE AMOUNTS ACTUALLY RECEIVED BY JAPET MEDICAL DEVICES®. JAPET MEDICAL DEVICES® SHALL NOT BE LIABLE FOR ANY INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES, LOSSES OR EXPENSES (INCLUDING BUT NOT LIMITED

TO LOST PROFITS) ARISING DIRECTLY OR INDIRECTLY FROM THE SALE, INABILITY TO SELL, USE OR INABILITY TO USE ANY PRODUCT. UNLESS OTHERWISE STATED IN THIS WARRANTY, ALL PRODUCTS ARE PROVIDED «AS IS» WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESS OR IMPLIED.

Registered trademarks

JAPET MEDICAL DEVICES® is a registered trademark of JAPET MEDICAL DEVICES®

Patents

The JAPET.W device is covered by several patents in France and internationally.

10. Other information

Symbols on the label and packaging :

 The symbol consists of the letters "MD" in a bold, sans-serif font, enclosed within a square border.	Indicate that the article is a medical device.
 A warning symbol consisting of an exclamation mark inside a triangle.	Symbole indicating potentially dangerous situations and actions.
1. 2. 3.	Detail of tasks for the correct use of the product.
 A silhouette of a factory with three chimneys.	The name and address of the manufacturer.
 A stylized outline of a factory with three chimneys.	The date of manufacture of the device.
 The letters "SN" in a bold, sans-serif font, enclosed within a square border.	Serial number of the device.
 The letters "REF" in a bold, sans-serif font, enclosed within a square border.	The device manufacturer's reference.

Consult the user guide carefully

CE marking of the device

DC power supply

10,8V

Voltage of the power supply supplied by the battery

3A



Power supply current delivered by the battery

Applied part type BF.

Symbol indicating the low and high temperature limit for storage and transport.

	Fragile, handle with care.
	Protect from moisture.
	Moisture restriction.
	Safe from the light
IP21	Protection class
	Separate system for waste collection
	Do not wash the device
	Batteries
	Charger
	Cables

Version history of the user guide

Version	Date	Modifications
1.0	10/2020	Creation based on the Atlas user guide.
2.0	11/2020	Historical and counter-indication corrections
3.0	02/2021	Translation into German and addition of a contraindication. Updated according to regulatory references.
4.0	03/2022	Change of address, regulatory references, belt tightening device, technical specifications, service contact.
5.0	07/2022	Clarification of the indication, addition of adverse effects, vigilance, clinical performance, modification / addition of symbols

Contact Informations

Japet Medical Devices®
F3A
12, avenue Pierre Mauroy
59120 LOOS
FRANCE

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